

Product datasheet for **VC102295**

ORF005L (NC_003494) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ORF005L (NC_003494) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ORF005L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC102295 represents NCBI reference of NP_612227 with codon optimized for human cell expression Red =Cloning site Blue =ORF Green =Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGTATGTGCCCTGTTGGAATCACGCCTTGAGAACCACAATGAGCTGCGTAGTTTCTACCCCTTGTC
AGACAATTATCTCCGTGTCTCGAGCCGCTATGTGTGCTGTGCCTTTGCGCCCCCTTCTCATCTACGGG
TGCTACCCGGCTTAGTTCCCATGAGACTATCGCACTGGCCATCGCTTCTAAAACCATGTTGACTGTTCTG
TTGGACCTGGACCACACGCTTATCGAATCCGTTAGTACCGGTAAGCTGAAGACAGCTAACGTGCAGAGAG
GACCGGGAAATGTTCCGACCCATGGAGCGTCCGCACTGCGACTTGCTAACGATATCGCTCGCTTCAATG
GCACCTCATCCCTGGTTACGTGGTCTGTTTGGCGCCTCACGCCAGAGCTTTTATAGCCCTCTTGTTGCA
GCCGGCTACAGGGTGGCGTCTGGACCGCCGCTCTGCCGTTATGCACGAGATATCCTGGCTCTGCTCC
ACCTGGACCGAATGCCCTTGCTACTTCATGTATGATGTGCACGACCAGACGGGTGTAAGCCACTGCA
TACCCTCTCGAGTCAAACGACATGGTAGTGATCGTCGACGACAACGATGCTGTGAAGGCCACACAGCCT
CACAGATGCGTGAAGTCCGCCCTTCTATGTTTCCACCCTGAGCAGATACCGCCCTCATGGACGCGAG
CTAGTGTGATCATAAACGATTTCCAGCCTTCCAGCTCCAGCAGCAGTACAAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC102295 representing NP_612227
 Red=Cloning sites Green=Tags

MRYVPCWNHALRTTMSCVVSTPCQTIISVSRAAMCRCVLCAPFSSTGATRLSSHETIALAIASKTMLTVL
 LDLDHTLIESVSTGKCLKTANVQRGPGNVPTHGASAVRLANDIARFEWHLIPGYVCLRPHARAFIALLFA
 AGYRVGVWTAASAGYARDILALLHLDRMPLSYFMYDVHDQTGGKPLHTLFE S NDMVVIVDDNDAVKATQP
 HRCVQVPPFYVSTPEHDTALMDAASVIINDFQPSQLQQQYK

TRTRRLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_003494

ORF Size: 753 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NC_003494.1, NP_612227
RefSeq ORF: 753 bp
Locus ID: 935358
MW: 27.6 kDa