

Product datasheet for **VC102304**

ORF014R (NC_003494) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ORF014R (NC_003494) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ORF014R
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC102304 represents NCBI reference of NP_612236 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATATTTATTACCAAGATGTATTCCTCCTGGAGATCGGAGAGGTGCTGGCTATTGGCAGGCCGACA
ACACCGTGTCTTTTATCCATGATCGCGACAATGAAGTCTACGTCGTGGTCCAGATTTCCGGGGCTGCGT
CCGCGTTATCGCCATGGCCGGACAAGACACCCACGCCGGCAAGCTGCGCGCTGTCTTTACACTGCGGCA
CACGCATACGGCCCCAACTTACTATGCCTCTGGAGTTTGCCAGACCATATATGGCCACAATGGTTTTG
TCAACTTCCATAGACCTGGGGAAGTAGGCATGTGCCCTGCTTATATGTCCGCACAACACTCCGCTTTGGT
TGGCCATACAGCCTACCACTCTGAGTGGAGCAACGCTGATACAGGTCAACGCTGTTATGAGCCTTCTG
CAGGTATCCCACTTCAAGGAAATGGTTGCCTATGACATCGTGGGGCCGGACCCCTCCATGTGTTGACCG
TTGAAGATGACCACATCTACCTGGTTAGTTGTATGTCCATGCCAGAATATGCAGATGCCAGCCCGCTGA
TACCAGATCCGCTTTGAGTACGCGGTGCGCATCGACATCCGCGACAAAATAGTGAGGTGGCAGGATGTT
TCAGTTCTGGTGGTCCCTCACCGCACCGGTCTGAGTGTATAGTGACAAGCGCTCAAGGCCGCCCGAGC
CTATCCTGATCCCATGGGATACTATCCAGGGGGAGTGTCCGTGACTACCCATGTTAAAGCTCCCCACAC
TCAGGGACTGATGCATGTGACATGTATAATATTAGCTCTCAACTGGCCACAATACCAAAGGTCAGGAT
ATTTTCGAAACTATGCCTCTCACACTTCGGACCGAGACGTTTCTCTGGACTGATAGCGAGTTTCGATA
CTGACTCCGAGGTGGCGATCCACCCGGCCACAAGTACAGTGCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC102304 representing NP_612236
Red=Cloning sites Green=Tags

MIFITKMYSLLEIGEVLAIGRPDNTVSFIHDRDNEVYVVVQISGGCVRVIAMAGQDTHAGKLRACLTYAA
 HAYGPKLTMPLEFARPYMAHNGFVNFHRPGEVGMCPAYMSAQHSALVGHTAYHTLSGATLIQVNVMSLL
 QVSHFKEMVAYDIVGAGPLHVLTVEDDHIYLVSCMSMPEYADASPPDTRSAFEYAVGIDIRDKIVRWQDV
 SVLVVPHRTGLSVIVTSAQGRPEPILIPWDTIQGGVSVTTHVKAPHTQGLMHVTCNLISSQLATITKGQD
 IFETMPLTLRDRDVS�DTSEFDTDSEVGDPPAAQSTVL

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NC_003494

ORF Size: 957 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_003494.1](#), [NP_612236](#)
RefSeq ORF: 957 bp
Locus ID: 935367
MW: 34.8 kDa