

Product datasheet for **VC102329**

ORF039R (NC_003494) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ORF039R (NC_003494) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ORF039R
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC102329 represents NCBI reference of NP_612261 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACCGGTTTACCCTTACCAGACTGCAGCTGGCCTACATGCTTCTCAAGCGCAACATGCAACCCGTCC
CCCAAACCTTCGATGCTAACGCGGATCCCCCGAACTGCTGCAAGCTATGTATAACCAGCCATGTGATAG
TCTCGGAGAGGTTCCGAGACAACACTACTGCGCCTGTGCAACATACTGACTCGATGAACAGCTCTTGGCT
GCGACAGATGACTGACAGACAGTATCTGTGGCTCTGCGCCATACACATATGCCTTACCGCCCTCCCC
GAGTGGATGTGTAGCTGCCGCGGTGGTCTGGCGACACAGCGCCTCATCGACAAGATGCGCGGATTTGA
ACAAATGGCAGCTTCTGTCTCAGAGCAGCTTCTGACCCATCAGCACTCAGATTCTCTTAAGCCTCTCTAC
GAGCTGCTGGATACTGTGCTGCCAAACTCTCCACCGTCAACAGGAACTTGGCCGAATCTGACGCACTGC
TGGATGCCGAATACGTGCTCAAAGATACGCTGGATGCCGACTTACCGAGACTATCCACTTGACTTATAA
GCCCCACTGGTGGTGGCCCTGGCTTTTGCCAAGCAGACCGGTGCAACCGCGCCTTTGAGAGTGCACCA
ATTCAGTCCCCAAGACCACAGCATGTTGGCCGTTTATGGCACTGTAATTTGTGTAGGGTCACTTGACG
TAGCATTGGTCTCGGCCAACCTTTGCTACATGCTGCGGACAACAATGAATGGCGTCGGCGCAGTGCC
CACCGCGGCTATAGTACACCACAAACAACAATCAGAATTGTGCTGATCGCTCCCAACCATTTACCCAG
GTCAACTTGACGGCTGTGATCCCTTTTAGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC102329 representing NP_612261
Red=Cloning sites Green=Tags

MDRFTLTRLQLAYMLLKRNMQVPVQTFDANADPPELLQAMYNQPCDSLGEVRRQLTAPVEHTVLDEQLLA
 ATDDTDRQYLWLCATISHMPYRPPRVSDVSAAVVLATQRLIDKMRGFEQMAASVSEQLLTHQHSDSLKPLY
 ELLDVTVAKLSTVNRNLAESDALLDAEYVLKDTLDARLTETIHLTYKPLVPLAFKQGTGRTGAFESAP
 IQVPKTHSMLAVYGTVNCVGSLDVALVSGPTFRYMLRRTMNGVGAVRTGGYSTPQTTIRIVLIAPNHLTQ
 VNLTAVIPFR

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_003494

ORF Size: 870 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NC_003494.1, NP_612261
RefSeq ORF: 870 bp
Locus ID: 935355
MW: 32.1 kDa