

## Product datasheet for **VC102368**

### ORF078R (NC\_003494) Virus Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ORF078R (NC\_003494) Virus Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** ORF078R  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >The Viral ORF clone VC102368 represents NCBI reference of NP\_612300 with codon optimized for human cell expression  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGATATCCCAGCTGCGGATGCACACTTTCACCCCACTTTCAGTGAGGAGAACCTTAGGGTTGACGGGG  
ATACTTTCGCAACCGTCATGTTGGAGCCCCCGACCACAGTACAATAAACATTGGACCTTCTTCACTAC  
CGACCATCTGCTTTTGTGCTATACCTGGTGCTTGAACAAGCGACGTTAACACCAACACAATATGATC  
GTGATGGCATCAGAGAACGGACTGTCAAAGACCCATCACATAGTCCCTGCCGTGGTGCAGGATAAGGTGT  
CCCTGTTTAGTGGTTTACATGATGCGAGTTGCTAAAGACCAGGAATCCATGTATAAGGTGAATATCTT  
TATCACGAGCACCGCCGAGGCACTCTCAAATGCTCGGCCGCAATGTTACCATGGTCTTTCTAGGCTG  
CATCAGCTACCAAGAGAACTGCAAAGCCGAGCTCATTGCT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >VC102368 representing NP\_612300  
**Red**=Cloning sites **Green**=Tags  
 MISQLRMHTFTPTFSEENLRVDGDTFATVMLEPPDHSTINIGPFFTTDHLCCYTLVLETSVNTKHNMI  
 VMAENGLSKTHHIAPAVVQDKVSLFSGFHMMRVAKDQESMYKVNIFITSTARGTLKMLGGNVTMVL SRL  
 HQSTKRATAKELIA

**TR**TRRLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



**Cloning Scheme:**


**ACCN:** NC\_003494

**ORF Size:** 462 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NC\\_003494.1](#), [NP\\_612300](#)

**RefSeq ORF:** 462 bp

**Locus ID:** 935385

**MW:** 17.2 kDa