

Product datasheet for **VC102375**

ORF085R (NC_003494) Virus Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | ORF085R (NC_003494) Virus Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | ORF085R |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| ORF Nucleotide Sequence: | >The Viral ORF clone VC102375 represents NCBI reference of NP_612307 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s) |

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAACATAGCATTGTATAGTACTATGCGCTACCGAACAAGATCAGAAAGATGGACTGCATTTGCCACC
TTGATTGCTTGTTCACCAACCCCTTCACCAAGTGTCTATGTTCTTCCCATCAACACACCCCTTCATGCC
TCAAGCCAAGACTATGTGGTGTGACGATTACCTCACACGACGCCGGCTCGAGCTGATGCCCTCCACCT
CCGCCACCAACCAAGCCAGTCACTGAGGAGTCACTTTCAGGTGGCCAGCCTCCCTCCTCCTT
CACCCTCCCAAGGCCAGCGAGTGTACGGAAGAGAGCAGCAGCGCTGGACCTCTCAGGCAGTGGTTC
CAGTGCTGACCATGAATCCCGGCGCAAAAAGAGGCACCGAAGGCGGCATCCTTATTCCCGCTCAGAGCGA
ATTTGCCATAATAATCCGAGAGGGTGCAGCGGATGCAAGTCAAAGATGATATGGACTGCCTTGCCACGA
TCATCGAACAGTGCAGACACCCGTGCAACGACGCACAAAGAGAGAAATCATTGCGTCAGCAATTCAGT
TTTGCTTAGCATGCAACAGAAGCTGAGGCAAGCTGCTCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

| | |
|-------------------|--|
| Protein Sequence: | >VC102375 representing NP_612307 Red=Cloning sites Green=Tags |
|-------------------|--|

MNIALYSTMRYNKIRKMDICICHLDFNPFKSCSSHQHTPSCPQAKTMWCDDYLRRRLELMPPPP
PPPPKASQCTEESSGGQPPPPSPLPKASECTEESSGWTSSGSGSSADHESRRKRGTGRRHPYSRSE
ICHNNSERVRRMQFKDDMDCLATIEQCIGHTVQRRRTKREIIASAIQCLLSMQQLRQAAQ


TRRRLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_003494

ORF Size: 600 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_003494.1](#), [NP_612307](#)

RefSeq ORF: 600 bp

Locus ID: 935352

MW: 22.7 kDa