

## Product datasheet for **VC102384**

### ORF094L (NC\_003494) Virus Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ORF094L (NC\_003494) Virus Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** ORF094L  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >The Viral ORF clone VC102384 represents NCBI reference of NP\_612316 with codon optimized for human cell expression  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTATCTTAGGGTGGTCGTTTGTTC AATCATGCACTGTGCTGTACTCACATTCGGAGATGTTTGC GTGT  
CCATTAATAGTGACGGAGCGGTCCACGGCCCTGGGCTGCTGACATAGGAGCGTACGCAGTGGCCATGGC  
CGCAGCTCACACCCACTACCAGCTGCAGTCAAATAATCTCTGGTGCAGCAATTGTGCGCCGGTGGGGCC  
GCTGTAGGCCACGCATTTGAGTGTGCTTACGCCATAGAGTGTGAGTCTTTCGGGCGGTTGATATCATGT  
GTTTTTGCAGCATACACAACACACTGCCGTGACAATTCAGTGCAGCGGCAGACTGTGTGCCAAAGGCAC  
GCTGCAAATCCCTAGTATGACTAGGATCGTCAGAGTAGTTGCTGGAGGACACATCTCCACCAAGTGTGG  
TGGAACGCTATGGATAACACAGACGACGTTTCGGACTCTCGGAAGAGGAACACATATCTGCATGCCTTGCA  
TGCACATG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >VC102384 representing NP\_612316  
**Red**=Cloning sites **Green**=Tags

MYLRVVVCSIMHCAVLTFGDVCSINS DGAVHGPWAADIGAYAVAMAAAHTHYQLQSNLWCQQLCAGGA  
AVGHAFECAYAIECESFGRVDIMCFDDTQHTAVTIQCSGRLCAKGLQIPSMTRIVRVVAGGHI SPPVW  
WNAMDNTDDVRTLGRGTHICMPCMHM

**TR**TRRLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NC\_003494

**ORF Size:** 498 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NC\\_003494.1](#), [NP\\_612316](#)

**RefSeq ORF:** 498 bp

**Locus ID:** 935416

**MW:** 17.9 kDa