

## Product datasheet for **VC102425**

### NS1 (NC\_007364) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	NS1 (NC_007364) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NS1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC102425 represents NCBI reference of YP_308673 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGACTCTAACACCATTACATCCTTTTCAGGTCGATTGTTACCTGTGGCACATCAGAAAGCTGCTCTCAA  
TGCCGACATGTGTGATGCACCTTTTGACGATAGACTGCCGAGAGATCAGAAAGCCCTGAAAGGTAGGGG  
CAGCACACTCGGCCTTGATTTGAGGGTTGCTACAATGGAAGGCAAGAAGATCGTGGAGGATATCCTGAAA  
AGTGAAACAAATGAAAATCTGAAGATTGCTATGCCAGCAGTCCAGCACCTCGATACATCACCGACATGA  
GCATCGAAGAGATGAGTCGCGAGTGGTACATGCTGATGCCGAGACAGAAGATCACCGAGGGTTGATGGT  
GAAGATGGATCAGGCAATAATGGACAAGAGGATCATTCTGAAAGCAAATTTCTCCGTGTTGTTCCACCAG  
CTGGAGACACTCGTGTCCCTCCGCGCATTCACTGAATCCGGGGCTATCGTGGCCGAGATCTCCCCATAC  
CAAGTGTGCCTGGCCATTTACCAGAACGCTGAAAAACGCTATTGGGATACTGATTGGCGGACTGGAATG  
GAATGACAATTCAATAAGAGCTTCTGAAAATATCCAACGGTTCGCGTGGGGGATCCACGATGAGAACGGA  
GGGCCAAGCCTTCCACCCAAACAAAAGCGATACATGGCCAAAAGAGTGGAGAGCGAAGTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NC\\_007364.1](#), [YP\\_308673](#)

**RefSeq ORF:** 690 bp

**Locus ID:** 8656648

**MW:** 26.1 kDa