

## Product datasheet for **VC102449**

### HA2 (NC\_007374) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	HA2 (NC_007374) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	HA2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC102449 represents NCBI reference of YP_308874 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGGGCTCTTCGGAGCAATTGCAGGCTTCATTGAAGGAGGCTGGCAAGGTATGGTCGACGGCTGGTACG  
GATATCATCACTCTAATGACCAGGGATCTGGTTACGCGCGGACAAAGAAAGCACACAGAAAGCGTTCAA  
CGGGATCACTAATAAGGTCAACAGTGTGATCGAGAAGATGAATACTCAGTTTGAGGCAGTCGGAAGGAA  
TTCTCCAACCTGGAAAAAAGATTGGAGAATCTGAACAAAAAGATGGAAGATGGGTTTCTGGATGTGTGGA  
CCTACAATGCCGAAGTCTGGTGTGATGGAGAATGAGAGGACTCTGGACTTCCATGATTCTAACGTCAA  
AAACCTGTACGACAAGGTTTCGCATGCAGCTCAGGGATAACGTAAGGAGCTCGGGAATGGATGTTTTGAG  
TTTTACCATAAATGCGACAACGAGTGCATGGATTCTGTTAAAAACGGCACGTATGATTATCCTAAATATG  
AAGAGGAGAGCAAGCTGAACCGAATGAGATTAAGGGAGTGAAGCTGAGTTCAATGGCGTTTATCAAAT  
ACTCGCTATATACGCTACAGTTGCAGGGAGCCTGAGTCTGGCCATTATGATGGCGGCATATCCTCTGG  
ATGTGTTCTAATGGTTCCTGCAGTGCAGGATATGCATC

**ACGCGT**ACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >VC102449 representing YP\_308874  
 Red=Cloning sites Green=Tags

MGLFGAIAGFIEGGWQGMVDGWYGYHHSNDQSGSYAADKESTQKAFNGITNKVNSVIEKMNTQFEAVGKE  
 FSNLEKRLENLNKMKMEDGFLDVWTYNAELLVLMENERTLDFHDSNVKNLYDKVVMQLRDNVKELGNGCFE  
 FYHKCDNECMDSVKNGTYDYPKYEEESKLNRNIEKGVKLSMGMVYQILAIYATVAGSLSLAIMMAGISFW  
 MCSNGLSQCRICI

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NC\_007374

**ORF Size:** 669 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NC\_007374.1, YP\_308874  
RefSeq ORF: 669 bp  
MW: 25.3 kDa