

Product datasheet for **VC102478**

NP (NC_002019) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NP (NC_002019) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>The Viral ORF clone VC102478 represents NCBI reference of NP_040982 with codon optimized for human cell expression
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCAGCCAGGGCACCAACGGTCATACGAGCAGATGGAAACCGACGGCGAAAGGCAGAATGCGACCG
 AGATTAGGGCTTCTGTGGGCAAAATGATCGGCGGAATTGGACGGTTTTACATACAATGTGTACCGAGCT
 GAAGCTGTCTGATTATGAGGGTAGGCTCATTGAGAACAGCTTGACTATAGAACGAATGGTCTGTCTGCC
 TTTGATGAGAGAAGAAATAAGTACCTCGAAGAATCCTTCAGCAGGGAAGGACCCTAAGAAGACAGGGG
 GTCCAATTTACCGCCGCTCAACGGGAAGTGGATGAGAGAATCATACTGTACGACAAGGAAGAGATTAG
 AAGAATCTGGAGACAAGCCAACAATGGTACGATGCAACCGCCGGGCTCACCCACATGATGATCTGGCAT
 AGTAACCTCAACGATGCCACCTATCAGAGAACACGAGCTCTGGTACAGAACGGGGATGGACCAAGGATGT
 GCTCCTTGATGCAGGGAAGTACCCTCCCAAGGAGAAGCGGGCAGCTGGAGCCGAGTGAAGGGCGTAGG
 TACAATGGTTATGGAGTTGGTGCAGATGATAAAGCGGGGAATCAACGACAGGAACTTTTGGCGAGGTGAA
 AATGGGGCAAACTCGAATTGCGTACGAGAGGATGTGTAATATCCTCAAGGGCAAAATTTCAAACAGCTG
 CTCAAAAAGCGATGATGGATCAGGTACGCGAGAGTACGACCCCGGCAACGCCGAGTTTGAAGATCTCAC
 TTCTTGCACGGAGCGCCCTGATCCTGAGGGGATCTGTTGCCACAAAAGCTGCCTTCCCGCTTGGCTT
 TACGGGCCAGCCGTCGCGTCTGGGTATGATTTTGAGAGGGAAGGCTATTCCTGGTGGGCATCGACCCCT
 TCCGGTTGTTGCAGAACAGTCAGGTATACTCTCATCAGGCCTAACGAGAACCCCGCCACAAGTCCCA
 GCTTGTGTGGATGGCCTGCCACTCCGCAGCATTGAGGATCTCCGAGTGTGTCATTATCAAAGGCACG
 AAAGTCGTCCACGAGGTAAAGTGTCCACTCGGGGAGTGCAGATTGCAAGCAATGAAAACATGGAGACAA
 TGGAAAGTTCAACCTTGGAACTGAGGTCAGATACTGGGCAATCCGGACCAGGAGCGGGGGCAACACAAA
 CCAGCAGCGAGCTTCTGCCGGGCAAAATAAGCATACAGCCTACCTTCTCAGTACAGCGCAACCTCCCATT
 GACCGCACACAGTCATGGCCGATTTACTGGGAATACCGAAGGGAGGACCTCCGATATGAGGACCGAGA
 TAATTCGCATGATGAAAAGCGCTCGCCCTGAAGACGTGAGCTTTCAGGGCGGGGAGTGTTCGAGCTGTC
 CGATGAAAAGGCCGTTCTCCATTGTTCCAAGCTTGGACATGTCCAACGAAGGGTCATACTTTTTTGGT
 GACAATGCCGAGGAATATGATAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>VC102478 representing NP_040982
 Red=Cloning sites Green=Tags

MASQGTKRSYEQMETDGERQNAEIRASVGMIGGIGRFYIQMCTELKLSDYEGRLIQNSLTIERMVLSA
 FDERRNKYLEEHP SAGKDPKKTGGPIYRRVNGKWMRELILYDKEEIRRIWRQANNGDDATAGLTHMMIWH
 SNLNDATYQRTRALVRTGMDPRMCSLMQGSTLPRRSGAAGAAVKGVGTMVMELVRMIKRGINDRNFWRGE
 NGRKTRIAAYERMCNILKGFQTAQKAMMDQVRESRDPGNAEFEDLTFLLARSALILRGSVAHKSCLPACV
 YGPAVASGYDFEREGYSLVGIIDPFRLQNSQVYSLIRPNENPAHKSQLVWMACHSAAFEDLRLVLSFKGT
 KVVPRGKLSRTRVQIASNENMETMESSTLELRSRYWAIRTRSGGNTNQQRASAGQISIQPTFSVQRNLFP
 DRITVMAAFTGNTEGRTSDMRTEIIRMESARPEDVSFQGRGVFELSDEKAASPIVPSFDMSENGSYFFG
 DNAAEYDN

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/ja2337_h02.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:



ACCN: NC_002019

ORF Size: 1494 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_002019.1](#), [NP_040982](#)

RefSeq ORF: 1494 bp

Locus ID: 956531
MW: 56.2 kDa