

## Product datasheet for **VC102527**

### anchored capsid protein C (NC\_012532) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	anchored capsid protein C (NC_012532) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	anchored capsid protein C
Synonyms:	ZIKV_gp1, ZIKV
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC102527 represents NCBI reference of YP_009227206 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAAGAATCCTAAGGAGGAAATTAGACGCATCAGAATCGTAAATATGCTCAAGCGAGGAGTCGCCCGG  
TTAACCCACTGGGAGGCCTGAAACGACTGCCTGCAGGACTGCTCTGGGACATGGACCAATTAGAATGGT  
GTTGGCTATATTGGCCTTTCTGAGATTTACAGCCATCAAACCCTCTCTCGGACTGATTAACCGCTGGGGA  
AGTGTAGGAAAAAAGAAGCAATGGAAATTATAAAAAAGTTAAAAAAGATCTGGCTGCCATGCTGCGCA  
TAATTAACGCTAGAAAAGAAAAGACGACGGGTGCCGATACAAGTATAGGAATTATTGGACTTCTCTT  
GACAACCGCTATGGCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>VC102527 representing YP_009227206 Red=Cloning sites Green=Tags
-------------------	---

MKNPKEEIRRIRIVNMLKRGVARVNPLGGLKRLPAGLLLGHGPIRMVLAILAFLRFTAIPSLGLINRWG  
SVGKKEAMEIIKKFKKDLAAMLRIINARKERKRRGADTSIGIIGLLLTAMA


**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:	Sgfl-MluI
--------------------	-----------



[View online »](#)

Cloning Scheme:



ACCN: NC\_012532

ORF Size: 366 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC\\_012532.1](#), [YP\\_009227206](#)

RefSeq ORF: 366 bp

MW: 13.6 kDa