



Product datasheet for VC102575

SARS-CoV-2 Nsp5 Gene Tagged ORF Clone (Native Sequence)

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | SARS-CoV-2 Nsp5 Gene Tagged ORF Clone (Native Sequence) |
| Tag: | mRFP |
| Symbol: | NSP5 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-mRFP (PS100041) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| Restriction Sites: | SgfI-MluI |
| Cloning Scheme: | |

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

| | |
|-----------------|---|
| ACCN: | NC_045512 |
| ORF Size: | 966 bp |
| OTI Disclaimer: | The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag. |



[View online »](#)

| | |
|------------------------|---|
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NC_045512.2 , YP_009725301.1 |
| RefSeq ORF: | 966 bp |
| MW: | 35.6 kDa |