


Assay Solutions for Immunology

Ready- to-use | For Multiple Sample Types | Key Targets Available

Product Highlights

Human CD44 ELISA Kit



Sample Types:

Human serum, plasma, cell lysates, tissue homogenates and other biological fluids.

Sensitivity:

4 pg/ml

Range:

15.6 pg/ml - 1000 pg/ml

Catalog #.: EA200067

Human MAGEA3 ELISA Kit

Sample Types:

Human serum, plasma and other biological fluids.


Sensitivity:

12.9 pg/ml

Range:

62.5 pg/ml - 5000 pg/ml

Catalog #: EA200011



ELISA Kits & Matched ELISA Antibody Pairs

Over 2,900 ready-to-use ELISA kits for the detection of key biomarkers in immunology, neurology, viral research areas.

ELISA Kit Targets	Catalog Number	Sensitivity
NEW Human CD4	EA200072	47 pg/mL
NEW Human CDC1	EA200073	84 pg/ml
NEW Human MAGEA3	EA200011	12.9pg/ml
NEW MAGEA4	EA200022	33.5 pg/ml
NEW AAV1 Capsid	EA200068	9x10 ⁶ capsids/ml
Human IL-6	EA200005	1.28 pg/ml
Human CD-80	EA200057	0.058 ng/ml
Human HSD17B13	EA200004	0.1 ng/ml
Human Huntingtin	EA200074	57 pg/ml

Over 160 matched ELISA antibody pairs, ready for a wide variety of assay development platforms.

ELISA Pair Targets	Capture Antibody	Detection Antibody
NEW Tau (pT217)	TA600577	TA700576
NEW Tau (pT181)	TA600576	TA700576
NEW Tau (total)	TA600579	TA700579
CD4	TA600277	TA700298
MAGEA4	TA600567	TA700567
TNF-alpha	TA600026	TA700026
IL4	TA600022	TA700022
CD19	TA600284	TA700286
BUB1B	TA600521	TA700531
CD33	TA600374	TA700374
IDH1	TA600198	TA700199
IL10	TA600017	TA700017
SERPINA1	TA600007	TA700007
VCAM1	TA600239	TA700239
TNFRSF1A	TA600021	TA700021

For full offerings on ELISA kits and validated antibody pairs, visit www.origene.com

OriGene Headquarters
9620 Medical Center Dr
Rockville, MD 20850
+1-301-340-3188
custsupport@origene.co
techsupport@origene.co
sales@origene.com

OriGene EU
SchillerstraBe 5
32052 Herford, Germany
+49-5221-346060-0
info-de@origene.com

www.origene.com

